




3-ply Mask for medical use, **BLUE**  
Medical imaging compatible  
CE-Marking (Type IIR) - EN14683  
Made in France 

V 11.04.2022



**WITHOUT NOSE BAR**  
Compatible with medical imaging

Bacterial  
Filtration



Perfectly suited to all  
types of  
medical imaging

Comfortable  
and breathable

## The Advantages<sup>+</sup>

- 100% French manufacturing
- **WITHOUT NOSE BAR:** compatible with medical imaging
- Electrostatic filter fleece
- Very comfortable – Excellent respiratory comfort
- Without latex, silicone, PVC
- Without Graphene
- Manufacturing process guaranteeing the cleanliness of our products from the production line up to the packaging process



Made in France 

Category	Non-sterile Class I Medical device Single-use surgical mask (type IIR)
Standard and legislation	<ul style="list-style-type: none"><li>• EU Medical Devices Directive 93/42/EEC</li><li>• Complies with standard EN 14683:2019+AC:2019</li></ul>
Description	<ul style="list-style-type: none"><li>• Mask for medical use, non-sterile, for single use.</li><li>• Model : 1000 series</li></ul> <b>Specifications :</b> <ul style="list-style-type: none"><li>• <b>WITHOUT NOSE BAR:</b> compatible with medical imaging</li><li>• Size M</li><li>• Dimensions: M. 170 x 95 mm</li><li>• Bacterial filtration efficiency: see performance table</li><li>• Mask weight (+/- 5%):</li><li>• Type IIR : <b>M</b> 2,90 gr</li></ul>
Composition	<ul style="list-style-type: none"><li>• Inner and outer layers: polypropylene nonwovens spunbond .</li><li>• Intermediate layer: nonwovens meltblown in polypropylene.</li><li>• Elastic: Polyamide / spandex mesh ou polypropylene / spandex mesh</li><li>• French manufacturing lines – Controlled raw materials and products</li></ul>
Materiovigilance	The materials likely to come into contact with the wearer's skin do not present a risk of irritation or harmful effects to health.
Storage and preservation	Mask should be stored in a clean, dry, cool and well-ventilated place. Humidity less than 80%, temperatures between +5°C and +40°C. Maximum storage period: 5 years in their original packaging.



# 3-ply Mask for medical use, **BLUE**




## Medical imaging compatible

## CE-Marking (Type IIR) - EN14683

## Made in France

V 11.04.2022

Performances	EN 14683:2019	Bacterial filtration effectiveness (BFE)	Differential Pressure (Pa/cm²)	Resistance to projections (kPa)	Cleanliness microbial (cfu/g)
	Type				
	Type II-R	≥ 98 %	< 60	≥ 16	≤ 30
Duration of use	4 hours maximum (in one shift – continuous)				
Sector of use	<ul style="list-style-type: none"> <li>· General public</li> <li>· Professional, sectors for example : <ul style="list-style-type: none"> <li>· Medical sector not exposed to projections and splashes</li> <li>· Medical sector exposed to projections and splashes</li> <li>· Service, community, visitors</li> <li>· Infrastructure, construction</li> <li>· Light industry</li> <li>· Logistics transport</li> <li>· Medico-Social</li> <li>· Retirement home</li> <li>· Schools</li> <li>· Various exposed areas protecting the wearer against the health risk</li> </ul> </li> </ul>				
Packaging	Box of 50 units 30 boxes in one cardboard Packaged without bags Or 5 bags of 10 masks				



KLAP France SAS

Avenue des Châtelets – FR-24440 Ploufragan – Tel : +33 2 90 52 04 33 – hello@klapsolutions.fr – www.klapsolutions.com – member of m3 GROUPE